


# CELESTA light engine®

The Next Generation Comes To Light



lumencor®



## Integrated Array of Seven Solid-State Laser Light Sources

Confocal • TIRF • Super-Resolution • Optogenetics • FRAP • SPIM

The next generation of solid-state illumination is here. In Lumencor's CELESTA light engine, seven individually addressable solid-state laser light sources join forces with advanced electronic control systems to deliver unprecedented output power and performance.

The CELESTA light engine delivers approximately one watt of output power from each laser line (approximately seven watts total). Long-term power stability is sustained by active stabilization. An onboard photodiode continuously monitors the light output and generates a reference signal that is applied to the constituent sources in a feedback loop to maintain constant light output over time.

The laser outputs are refined by bandpass filters and merged into a common optical train directed to the light output port on the front panel. The light output port has a built-in adapter for connection to microscopes and other bioanalytical instruments through a SMA-terminated optical fiber.

The CELESTA features an advanced control system based around an onboard computer with an embedded

command library. This allows control using simple and intuitive text string commands sent to the light engine via USB/RS-232 or TCP serial protocols. These commands give access not only to the basic control functions of light source selection, on/off switching and output intensity adjustment, but also to an extensive panel of operating status reports and preference settings. A GUI resident on the onboard computer and viewed using a web browser via a LAN connection provides convenient access to many of the command library functions. CELESTA controls are also implemented in several image acquisition software packages. TTL trigger inputs are provided for all 7 output lines for applications requiring fast (10 microseconds) switching.

---

For more information on the CELESTA light engine, please contact us at [info@lumencor.com](mailto:info@lumencor.com). To receive a purchase quotation for a CELESTA light engine, please submit our [online quotation request form](#).

# CELESTA light engine®

The Next Generation Comes To Light



408

445


473

518

545

635

750



**NOTE:** 10X image of fern sorus autofluorescence  
Courtesy of Michael Davis, Nikon Instruments

## Features and Operating Characteristics:

Features	Details
Sources	7 solid-state laser sources
Wavelengths	Nominal center wavelengths 408, 445, 473, 518, 545, 635, 750 ± 2 nm
Bandpass Filters	Integrally installed bandpass filters for spectral output refinement
Output Power	~1 watt per laser line with optional optical power stabilization function
Light Delivery	SMA-terminated fiber
Safety Interlocks	Laser output contingent on manual (key) and remote (electronic) interlocks
Operational Control	Onboard computer with server/client architecture and embedded command library
Control Interfaces	Source selection, light output on/off and intensity via serial interface (RS-232/USB or TCP). Source selection and light output on/off via TTL
Software	Onboard GUI or PC-based image acquisition software
Power Requirements	220 W (24V DC/9.2A) power supply included
Warranty	36 months
Dimensions (W x L x H)	145 mm x 340 mm x 203 mm (5.7 in x 13.4 in x 8.0 in)
Weight	8.7 kg /19.1 lbs
Optional Accessories	7-channel breakout cable for TTL triggering. Light engine control pod[1]

[1] Control pod connects to light engine USB port and controls source selection, light output on/off and intensity settings.



## GET IN TOUCH

Lumencor, Inc.  
14940 NW Greenbrier Parkway, Beaverton, OR 97006 USA • T 503.213.4269 • [www.lumencor.com](http://www.lumencor.com)  
©2018 Lumencor, Inc. • Effective Date: 03/2018 • 54-10043