



light engines for a

BRIGHTER. GREENER. PLANET.

LIGHT **BYTES:** May 2016

Better Red Than Dead...

In addition to savings of \$500–\$1,000 per year in operating and maintenance costs, another major benefit of Lumencor's **SOLA light engines**® over metal halide lamps is their superior output in the red region (>620 nm) of the visible spectrum. Images of *selaginella* autofluorescence, excited from 617–645 nm and detected under identical conditions, illustrate the improved image data quality delivered by the **SOLA light engine**. Excitation at >620 nm provides advantages in terms of greater sample transparency, less generation of phototoxic reactive oxygen species and lower background signals compared to excitation at shorter wavelengths. Many red-excited fluorophores including Cy5, Alexa Fluor 647, ATTO 647N TO-PRO-3 and DRAQ5 have been developed to allow exploitation of these advantages.

Give your microscope a **SOLA** upgrade today and look forward to brighter days ahead!




FIGURE A
• SOLA SE Light Engine
• 50% maximum output




FIGURE B
• 120 W Metal Halide
• 100% maximum output

TECHNICAL DETAILS

- Single band exciter 631/28 (CWL/FWHM). Semrock FF660-DiO2 dichroic and FF01-692/40 emitter
- Nikon Ti microscope, 10X, 0.5 NA objective, Andor Zyla 5.5 camera, 50 ms exposures.
- Power output at sample plane (10X, 0.5 NA objective) SOLA SE (100%) = 94 mW. 120 W metal halide (100 %) = 27 mW.
- Sample: Autofluorescent longitudinal section of a strobilus from selaginella (heterosporous clubmoss).



Be sure to visit our site lumencor.com and [blog](#) for more information about Lumencor products and applications.

