

LIGHT **BYTES:** May 2018


New SOLA SE® Light Engines With Added Cy7 Excitation

In response to increasing demand for light engines with near-infrared (nIR) output, we are pleased to offer two NEW SOLA SE nIR light engine models:

SOLA SE U-nIR light engine (Part Number: 90-10432) provides the same features as the SOLA SE 365 model with spectral output extended to about 760 nm.

SOLA SE V-nIR light engine (Part Number: 90-10444) provides the same features as the SOLA SE 390 model with spectral output extended to about 760 nm.

For an introductory period, we are offering the SOLA SE nIR light engines **at the same list price** as the original SOLA SE 390 and SOLA SE 365.



Both SOLA SE nIR light engine models provide:

- ~ 4.0 W total light output through standard 3 mm diameter liquid light guide
- Light on/off and graduated intensity control via control pod or PC
- TTL electronic on/off
- No routine maintenance, no consumables, mercury-free

To receive a purchase quotation for a SOLA SE nIR light engine, please submit our [online quotation request form](#).

